


Supplement D

Example of anal pressure recording during continued low-frequency stimulation after the cessation of high-frequency stimulation



An example of a recording of anal pressure during unilateral proximal low-frequency stimulation (square wave, 20 Hz, 3 V amplitude, 20% duty cycle, no interphase delay) and distal high-frequency stimulation (square wave, 5 kHz, 10 V amplitude, 50% duty cycle, no interphase delay). After two seconds of proximal low-frequency stimulation, distal high-frequency stimulation was applied for ten seconds after which proximal low-frequency stimulation continued for two seconds.